



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Breast Milk Fridges

Breast Milk Fridge RLBM0219

This Labcold Breast Milk refrigerator has been specifically designed to store human milk in optimum conditions and features minimum/maximum temperature recording to assist with record keeping in accordance with NICE guidelines. The controller also automatically records temperature and refrigerator data. This data can be extracted by the use of a specially coded Micro SD (optional extra) and downloaded to a computer as a csv file.

With three adjustable shelves and an interior fan for even temperature distribution, this fridge is fitted with a lock and has audio and visual alarms to warn of temperature deviation and if the door is left open. For extra security, your Human Milk fridge can be fitted with a factory fitted digital lock which solves the problem of generic and lost keys.

Specifications

Product Code	RLBM0219
Gross Capacity	66 litres
Style	Bench top/wall mounted
UK Warranty	1 Year
External Dimensions (mm) HxWxD	735 x 450 x 510
Internal Dimensions (mm) HxWxD	588 x 331 x 226
Temperature Range	0°C to 5°C
Fridge Set Point	3°C

Specifications (cont.)

Internal Material	Moulded Plastic
External Material	Powder Coated Steel
Shelves	3
Net Weight	35kg
Door Style	Solid
Lock Type	Key lock
Digital Temperature Display	Yes
Digital Temperature Control	Yes
High/Low Temperature Alarm	Yes
Min/Max Temperature Recording	Yes
Integral Temp Data Logging	Yes
Automatic Temperature & Fridge Operation Data	Yes
Number of Temperature Probes	2
Internal Light	Yes
Self Closing Door	Yes
Door Open Alarm	Yes
Mains Failure Alarm	Yes
Alarm Battery Backup	Yes
Auto Defrost	
Fan Assisted Cooling	Yes
Remote Alarm Contacts	Yes
Access Porthole	Yes
Access Porthole Size	12mm
Eco Friendly Refrigerant	Yes
Power Consumption Kwh/24hrs	
Ambient Operating Temperature	16°C to 32°C

Features

- Fitted with a Touch Screen IntelliCold® controller
- 24/7 automatic temperature data logging
- Power failure alarm
- Low temperature alarm
- Off cycle auto defrost
- Second temperature probe
- Ammonia free
- Real time temperature display
- Optional micro SD card for transferring temperature data to your computer
- Can be wall mounted on integral brackets
- Controller and alarm battery back up and door open alarm
- High temperature alarm
- Minimum/maximum temperature recording and display
- Fan assisted cooling
- Environmentally friendly refrigerant
- Door lock with two keys
- Internal light
- FREE 1 year parts and labour warranty

Models

	Model	Gross Capacity	Door Style	Lock Type
	Breast Milk Fridge RLBM0219	66 litres	Solid	Key lock
	Breast Milk Fridge RLBM0224	66 litres	Solid	Key lock

Accessories



Micro SD Card



Data Logger
RMBL2019-M1



Digital Thermometer
and Alarm
RLAA5003-1