



LABCOLD™

medical & scientific refrigeration



RLDG0510 Glass Door Pharmacy & Vaccine Refrigerator

Powered by the unique IntelliCold™ Controller, this pharmacy refrigerator has superior temperature uniformity and energy consumption, plus extra functions designed to store temperature sensitive pharmaceuticals in optimum conditions and assist busy professionals with regulatory requirements.



The simple to use four button interface is clearly labelled and allows minimum/maximum temperatures to be displayed at the touch of a button. It also shows clearly which alarm is active and why, making remedial action easy. A clear LED display shows the internal temperature of the fridge to one decimal place.

For objective and accurate record keeping, once the warranty is registered, a complimentary SD card will be sent which can be used to download all the temperature data stored by the controller. This can also be used to extract short term data if, for example, there is a suspected temperature incursion; data can be downloaded and used to prove the viability of the contents. Extra reassurance is provided by controller battery back up which can record and alarm in the event of a mains power failure.



The external light switch on the IntelliCold™ Controller means that the light can operate independently of the door opening, meaning you can illuminate the contents when the door is closed.

We are confident that this pharmacy refrigerator will provide you with many years of fault free service. However, for extra peace of mind your RLDG0510 comes complete with a FREE 5 year parts and labour warranty!*

**CLAIM
YOUR
FREE** **5** **YEAR
WARRANTY**

* Terms and conditions apply, please visit the web site

Part Number	RLDG0510
Capacity (Litres)	150
Style	Under bench
Exterior Dimensions (HxWxD)	820 x 600 x 610mm
Interior Dimensions (HxWxD)	690 x 510 x 410mm
Energy Consumption	0.66 kWh/24h
Operating Temperature	+5°C
Interior Construction	Painted Corrosion proof Aluminium
Exterior Construction	Powder Coated Steel
Shelves	3
Internal Drawers Available	✓
Door Hinging	Right hand only
Door Style	Glass
Door Lock	✓
Access Port	✓
Digital Temperature Control	✓
Digital Temperature Display	✓
Internal Light	✓
High/Low Temperature alarm	✓
Minimum/Maximum Temperature Recording	✓
Up to 13 months temperature data storage	✓
Door Open alarm	✓
Mains Failure alarm	✓
Alarm Battery Back-Up	✓
Off cycle auto defrost	✓
Fan Assisted Cooling	✓
Remote Alarm Contacts	✓
Eco Friendly Refrigerant	✓
Ammonia Free	✓
DH Green Book, RPS/ GPhC, NaTHNaC and RCVS Guidelines compliant	✓

01256 705 570

www.labcold.com

