



# LABCOLD™

medical & scientific refrigeration



## RLDF18041 Solid Door Pharmacy & Vaccine Refrigerator

Part Number	RLDF18041
Capacity (Litres)	505
Style	Freestanding
Exterior Dimensions (HxWxD)	1710 x 765 x 725mm
Energy Consumption	1.8 kWh/24h
Operating Temperature	+5°C
Temperature Set Point Range	+2°C to +8°C
Interior Construction	Powder Coated Steel
Exterior Construction	Powder Coated Steel
Shelves	5
Door Hinging	Right
Door Style	Solid
Internal Light	✗
Door Locks	Standard
Access Port	Optional Extra
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Door Open Alarm	✓
Maximum Temperature Display	✓
Minimum Temperature Display	✓
Auto Defrost	✓
Fan Assisted Cooling	✓
Environmentally Friendly Insulation	✓
Eco Friendly Refrigerant	✓
Ammonia Free	✓

Robustly constructed, this **Labcold** refrigerator is designed to provide close temperature control for temperature sensitive medicines.

This 505 litre, large volume Pharmacy Refrigerator is ideal for busy pharmacies and wholesale pharmacists. It comes complete with all the technology required to ensure expensive stock is kept in optimum condition, including features such as integral temperature display and alarm and forced air cooling as advised by Public Health England in Protocol for Ordering, Storing and Handling Vaccines.



Digital Temperature Display



The digital display provides visual information on temperature display, with minimum and maximum temperature display. An audible and visual high and low alarm will sound if temperatures deviate from the correct levels, plus a door open alarm reminds users to close the door after use.

Completely compliant with all the regulations specified by both the Department of Health and the RPS/GPhC this refrigerator is the professional choice for safe temperature controlled volume vaccine and pharmaceutical storage.

01256 705 570

[www.labcold.com](http://www.labcold.com)

