



# LABCOLD™

medical & scientific refrigeration



## RLCF1520 Sparkfree Chest Freezer

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

**Warning: Do not store open containers of volatile substances in this freezer.**

**Did you know that 25% of chemical accidents are caused by improper storage?**

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security. This model comes with 2 storage baskets to assist with organisation and storage. Please ask at time of order for additional baskets. This chest freezer has manual defrost for greater operator control.

Part Number	RLCF1520
Capacity (Litres)	447
Style	Chest Freezer
Exterior Dimensions (HxWxD)	860 x 1300 x 730mm (Including Handle)
Interior Dimensions (HxWxD)	720 x 1180 x 530mm
Weight	68kg
Energy Consumption	1.52kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Aluminium
Exterior Construction	White Enamel
Lid Style	Solid
Door Lock	✓
Baskets	2
Basket Dimension (HxWxD)	200 x 250 x 500mm
Feet	4 (adjustable)
High Temperature Alarm	Visual
Manual Defrost	✓
Eco Friendly Refrigerant	✓

01256 705 570

[www.labcold.com](http://www.labcold.com)

