



## Labcold Cube Ice Makers



LITV - NDP20



**Labcold** Cubed Ice Makers produce cylindrical cubed ice that is both hygienic and clear so ideal for medical and laboratory applications. The unique cylindrical shape makes them slower to melt than normal cubes.

Robustly constructed from stainless steel with well-insulated plastic collection bins and a lime scale reducing spray system.

**Labcold** Ice Makers require both a water and electricity supply and a drain at least 150mm below the ice maker.

Installation of an anti-scale filter is advised in areas with medium or high limescale content (hard water).



Ice Dimensions  
38mm (h) x 35mm (w)



LITV-DP80

| Part Number                 | LITV-NDP20        | LITV-DP35H                               | LITV-DP80                                  | LITV-DP150                                  |
|-----------------------------|-------------------|--|--|---|
| Ice Storage                 | 6kg               | 20kg                                     | 35kg                                       | 75kg  |
| Daily Ice Production        | 23kg              | 43kg                                     | 107kg                                      | 150kg                                       |
| Style                       | Self Contained    | Self Contained                           | Self Contained                             | Self Contained                              |
| Exterior Dimensions (HxWxD) | 593 x 350 x 475mm | 894/944 (789 without legs) x 465 x 595mm | 1044/1199 (936 without legs) x 715 x 595mm | 1149/1199 (1044 without legs) x 860 x 700mm |
| Weight                      | 55kg              | 45kg                                     | 70kg                                       | 95kg  |
| kW/h                        | 28.2              | 20.8                                     | 17.3                                       | 28.6  |
| Water use - Litres per hour | 3.41              | 3  | 5.9  | 19.2  |
| Refrigerant                 | R290              | R290                                     | R290                                       | R290  |
| Exterior Construction       | Stainless Steel   | Stainless Steel                          | Stainless Steel                            | Stainless Steel                             |
| Automatic Production        | ✓                 | ✓  | ✓  | ✓   |
| Feet                        | N/A               | 4  | 4  | 4   |