



## Labcold Cube Ice Makers



LITV - NDP20

**Labcold** Cubed Ice Makers produce cylindrical cubed ice that is both hygienic and clear so ideal for medical and laboratory applications. The unique cylindrical shape makes them slower to melt than normal cubes.

Robustly constructed from stainless steel with well-insulated plastic collection bins and a lime scale reducing spray system.

**Labcold** Ice Makers require both a water and electricity supply and a drain at least 150mm below the ice maker.

Installation of an anti-scale filter is advised in areas with medium or high limescale content (hard water).



Ice Dimensions  
38mm (h) x 35mm (w)



LITV-DP80

Part Number	LITV-NDP20	LITV-DP35H	LITV-DP80	LITV-DP150
Ice Storage	6kg	20kg	35kg	75kg
Daily Ice Production	23kg	43kg	107kg	150kg
Style	Self Contained	Self Contained	Self Contained	Self Contained
Exterior Dimensions (HxWxD)	593 x 350 x 475mm	894/944 (789 without legs) x 465 x 595mm	1044/1199 (936 without legs) x 715 x 595mm	1149/1199 (1044 without legs) x 860 x 700mm
Weight	55kg	45kg	70kg	95kg
kW/h	28.2	20.8	17.3	28.6
Water use - Litres per hour	3.41	3	5.9	19.2
Refrigerant	R290	R290	R290	R290
Exterior Construction	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
Automatic Production	✓	✓	✓	✓
Feet	N/A	4	4	4