



Labcold Cubed Ice Makers



LITV-NDP20/DP35H

Labcold Cubed Ice Makers produce cylindrical cubed ice that is both hygienic and clear so ideal for medical and laboratory applications. The unique cylindrical shape makes them slower to melt than normal cubes.

Robustly constructed from stainless steel with well-insulated plastic collection bins and a lime scale reducing spray system.

Labcold Ice Makers require both a water and electricity supply and a floor level drain.



LITV-DP80/DP150

Part Number	LITV-NDP20	LITV-DP35H	LITV-DP80	LITV-DP150
Ice Storage	6kg	20kg	35kg	75kg
Daily Ice Production	20kg per day	35kg per day	81kg per day	150kg per day
Style	Freestanding	Freestanding	Freestanding	Freestanding
Exterior Dimensions (HxWxD)	601 (595 without legs) x 350 x 475mm	900 (795 without legs) x 465 x 595mm	1050 (945 without legs) x 715 x 595mm	1155 (1050 without legs) x 860 x 700mm
Interior Construction	Plastic Bins	Plastic Bins	Plastic Bins	Plastic Bins
Exterior Construction	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
Automatic Production	✓	✓	✓	✓
Door Hinging	Bottom inward folding	Bottom inward folding	Bottom inward folding	Bottom inward folding
Door Style	Solid	Solid	Solid	Solid
Feet	4	4	4	4