



Labcold Flake Ice Makers



LITV-IQ50/IQ85C

Labcold Flake Ice Makers produce crushed ice ideal for filling packs and moulding to fit any item or object. The ice produced is both hygienic and clear so ideal for medical and laboratory applications.

Robustly constructed from stainless steel with well-insulated plastic collection bins and a lime scale reducing spray system.

Labcold Ice Makers require both a water and electricity supply and drainage 150mm below the level of the ice maker.



LITV-IQ135

Part Number	LITV-IQ50	LITV-IQ85C	LITV-IQ135	LITV-IQ200
Ice Storage	20kg	20kg	40kg	80kg
Daily Ice Production	66kg	90kg	156kg	234kg
Style	Self Contained	Self Contained	Modular	Modular
Exterior Dimensions (HxWxD)	855 (790 without legs) x 465 x 595mm	855 (790 without legs) x 465 x 595mm	1355 x 515 x 550mm	575 x 515 x 550mm Silo 870 x 555 x 685
Weight	58kg	58kg	75kg	Ice maker 52kg Silo 50kg
Exterior Construction	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
Automatic Production	✓	✓	✓	✓
Feet	4	4	N/A	N/A