



## RLCF0720 Sparkfree Chest Freezer

Part Number	RLCF0720
Capacity (Litres)	215
Style	Chest Freezer
Exterior Dimensions (HxWxD)	865 x 725 x 655mm (including Handle)
Interior Dimensions (HxWxD)	697 x 606 x 536mm
Weight	46kg
Energy Consumption	0.71 kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Aluminium
Exterior Construction	White Enamel
Lid Style	Solid
Door Lock	✓
Baskets	1
Basket Dimension (HxWxD)	200 x 250 x 500mm
Feet	4 (adjustable)
Mains Indicator Light	✓
High Temperature Alarm	Visual
Defrost drain	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

**Warning: Do not store open containers of volatile substances in this freezer.**

**Did you know that 25% of chemical accidents are caused by improper storage?**

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security. This model comes with a storage basket to assist with organisation and storage. Please ask at time of order for additional baskets. This chest freezer has manual defrost for greater operator control.