



# LABCOLD™

medical & scientific refrigeration



## RLCF0720 Sparkfree Chest Freezer

Part Number	RLCF0720
Capacity (Litres)	200
Style	Chest Freezer
Exterior Dimensions (HxWxD)	865 x 725 x 655mm (including Handle)
Interior Dimensions (HxWxD)	650 x 520 x 450mm
Weight	43kg
Energy Consumption	0.43kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Steel inner
Exterior Construction	White Galv Steel
Lid Style	Solid
Lid Lock	✓
Baskets	1
Feet	4
Mains Indicator Light	✓
High Temperature Alarm*	Visual
Defrost drain	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

**Warning: Do not store open containers of volatile substances in this freezer.**

**Did you know that 25% of chemical accidents are caused by improper storage?**

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a lid lock for optimum security. This model comes with a storage basket to assist with organisation. Please ask at time of order for additional baskets.

This chest freezer has manual defrost for greater operator control and is backed by a one year UK parts and labour warranty.

\* Please note this product may only have a mains indicator light

01256 705 570

[www.labcold.com](http://www.labcold.com)

