



LABCOLD™

medical & scientific refrigeration



RLCF1520 Sparkfree Chest Freezer

Part Number	RLCF1520
Capacity (Litres)	527
Style	Chest Freezer
Exterior Dimensions (HxWxD)	865 x 1505 x 655mm (Including Handle)
Interior Dimensions (HxWxD)	697 x 1386 x 536mm
Weight	67kg
Energy Consumption	1.3kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Aluminium
Exterior Construction	White Enamel
Lid Style	Solid
Door Lock	✓
Baskets	3
Basket Dimension (HxWxD)	200 x 250 x 500mm
Feet	4 (adjustable)
Mains Indicator Light	✓
Defrost Drain	✓
Manual Defrost	✓

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this freezer.

Did you know that 25% of chemical accidents are caused by improper storage?

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security. This model comes with 3 storage baskets to assist with organisation and storage. Please ask at time of order for additional baskets or racking.

This chest freezer has manual defrost for greater operator control.

01256 705 570

www.labcold.com

