



# LABCOLD™

medical & scientific refrigeration



## RLCF1520 Sparkfree Chest Freezer

Part Number	RLCF1520
Capacity (Litres)	416
Style	Chest Freezer
Exterior Dimensions (HxWxD)	865 x 1505 x 655mm (Including Handle)
Interior Dimensions (HxWxD)	650 x 1300 x 430mm
Weight	67kg
Energy Consumption	0.72kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Steel inner
Exterior Construction	White Galv Steel
Lid Style	Solid
Lid Lock	✓
Baskets	2
Feet	4
Mains Indicator Light	✓
Defrost Drain	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

**Warning: Do not store open containers of volatile substances in this freezer.**

**Did you know that 25% of chemical accidents are caused by improper storage?**

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security. This model comes with 2 storage baskets to assist with organisation. Please ask at time of order for additional baskets or racking.

This chest freezer has manual defrost for greater operator control and is backed by a one year UK parts and labour warranty.

01256 705 570

[www.labcold.com](http://www.labcold.com)

