



# LABCOLD™

medical & scientific refrigeration

## RLCF2120 Sparkfree Chest Freezer

Part Number	RLCF2120
Capacity (Litres)	607
Style	Chest Freezer
Exterior Dimensions (HxWxD)	865 x 1705 x 655mm (Including Handle)
Interior Dimensions (HxWxD)	697 x 1586 x 536mm
Weight	77kg
Energy Consumption	1.52kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Aluminium
Exterior Construction	White Enamel
Lid Style	Solid
Door Hinging	Back (viewed from front)
Door Lock	✓
Baskets	3
Basket Dimension (HxWxD)	200 x 250 x 500mm
Feet	4 (adjustable)
Mains Indicator Light	✓
Defrost Drain	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓



The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

**Warning: Do not store open containers of volatile substances in this freezer.**

**Did you know that 25% of chemical accidents are caused by improper storage?**

**Source: US Chemical Safety and Laboratory Investigation Board**

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security.

This model comes with 3 storage baskets to assist with organisation and storage. Please ask at time of order for additional baskets.

This chest freezer has manual defrost for greater operator control and a defrost drain to assist with getting rid of any defrost liquids quickly plus it comes with a 1 year UK parts and labour warranty.

01256 705 570

[www.labcold.com](http://www.labcold.com)

