



LABCOLD™

medical & scientific refrigeration



RLCF2120 Sparkfree Chest Freezer

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this freezer.

Did you know that 25% of chemical accidents are caused by improper storage?

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security.

This model comes with 7 storage baskets to assist with organisation and storage. Please ask at time of order for additional baskets.

This chest freezer has manual defrost for greater operator control and a defrost drain to assist with getting rid of any defrost liquids quickly plus it comes with a 1 year UK parts and labour warranty.

Part Number	RLCF2120
Capacity (Litres)	416
Style	Chest Freezer
Exterior Dimensions (HxWxD)	865 x 1705 x 655mm (Including Handle)
Interior Dimensions (HxWxD)	630 x 1500 x 450mm
Weight	77kg
Energy Consumption	0.80 kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Steel inner
Exterior Construction	WhiteGalv Steel
Lid Style	Solid
Lid Lock	✓
Baskets	7
Feet	4
Mains Indicator Light	✓
Defrost Drain	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓

01256 705 570

www.labcold.com

