



# LABCOLD™

medical & scientific refrigeration

## RLCF2120 Sparkfree Chest Freezer



The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

**Warning: Do not store open containers of volatile substances in this freezer.**

**Did you know that 25% of chemical accidents are caused by improper storage?**

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security.

This model comes with 2 storage baskets to assist with organisation. Please ask at time of order for additional baskets.

This chest freezer has manual defrost for greater operator control and a defrost drain to assist with getting rid of any defrost liquids quickly plus it comes with a 1 year UK parts and labour warranty.

Part Number	RLCF2120
Capacity (Litres)	556
Style	Chest Freezer
Exterior Dimensions (HxWxD)	865 x 1705 x 655mm (Including Handle)
Interior Dimensions (HxWxD)	630 x 1500 x 450mm
Weight	77kg
Energy Consumption	0.80 kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Steel inner
Exterior Construction	White Galv Steel
Lid Style	Solid
Lid Lock	✓
Baskets	2
Feet	4
Mains Indicator Light	✓
Defrost Drain	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓