



## RLFF13246 Sparkfree Fridge Freezer

Part Number	RLFF13246
Capacity (Litres)	Fridge 228
Capacity (Litres)	Freezer 134
Style	Freestanding
Exterior Dimensions (HxWxD) mm	2002 x 595 x 598 (plus 30mm for LK models)
Weight	80.4kg
Energy Consumption	1.36kWh/24hr
Operating Temperature	Fridge +2°C to +10°C Freezer -18°C to -23°C
Factory Set Point	Fridge +4°C Freezer -20°C
Interior Construction	Moulded Plastic
Exterior Construction	White Enamel
Door Style	Solid
Door Hinging	Right
Door Lock	Lock Model RLFF13246-LK
Shelves	Fridge 4 glass
Drawers	Freezer 4 plastic
Feet	4
Temperature Display	Optional
High Temperature Alarm	Optional
Low Temperature Alarm	Optional
Defrost	Fridge Auto Freezer Manual

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

**Warning: Do not store open containers of volatile substances in this fridge freezer.**

Did you know that 25% of chemical accidents are caused by improper storage?

Source: US Chemical Safety and Laboratory Investigation Board



Designed specifically for laboratory safety with a sparkfree interior and an optional factory fitted door lock for optimum security (ask for lock model RLFF13246-LK).

This fridge freezer provides a large storage capacity on a small floor area with two refrigeration systems so you can control the temperature of the fridge and freezer separately.