



## TYPICAL MAINTENANCE CHECKS

At **Labcold** we know there is no such thing as a 'typical' maintenance, but these are just some of the checks we can do as part of an annual preventative maintenance inspection:

- ❄️ Full mechanical and electrical inspection. Testing of the whole unit; including doors, gaskets, fans and wiring.
- ❄️ Test and record alarm system function and connection to remote alarm output.
- ❄️ Inspect and calibrate chart recorders and replace battery.
- ❄️ Check temperature sensors and simulated load bottle.
- ❄️ Check compressor functionality.
- ❄️ Check and calibrate digital display for both air and load temperatures where applicable.
- ❄️ Provide advice on best practice and update training of operators if necessary.
- ❄️ Issue a comprehensive maintenance report
- ❄️ Issue calibration certificate.

**Labcold** is the UK's leading manufacturer of professional medical and scientific refrigeration and can uniquely offer customers Original Equipment Manufacturers (OEM) preventative maintenance contracts. All **Labcold** blood storage equipment has medical device certification so must be maintained to the level of best practice as defined by the MHRA. As a result, the regulatory and safety benefit of a Labcold maintenance contract has never been more relevant.

At **Labcold** we take pride in providing our customers with a service which lasts the life cycle of their purchase. We care about quality and safety which is why we offer a cost effective range of annual maintenance contracts designed to suit your individual budgets and requirements. There are three **Labcold** maintenance contract options:

ANNUAL:	STANDARD	PREMIUM	PREMIUM+
Maintenance Visit	✓	✓	✓
Inspection/calibration check	✓	✓	✓
5% off spare parts*	✓	✗	✗
10% off spare parts*	✗	✓	✓
Full UK coverage	✓	✓	✓
Free calibration	✓	✓	✓
2 Free Fault Visits	✗	✓	✗
4 Free Fault Visits	✗	✗	✓
No Claims Discount on renewal	✗	✓	✓
Priority call out	✓	✓	✓

\*if fitted by a Labcold engineer

If your requirement is different we can adapt a contract to suit you, so why not contact us for a quote and more details?

We can also provide you with GMP compliant temperature mapping which can form part of IQ/OQ/PQ procedures. In addition, we have a trained service team who can answer your questions, free of charge by phone or e-mail.