




RVPF2134MD 505 litre Upright Plasma Freezer

The RVPF2134MD is the largest plasma freezer manufactured by Labcold. Like all Labcold plasma freezers, it has been designed to operate at a temperature that maintains plasma in optimum condition whilst minimising the risk of bag cracks appearing on thawing and is a certified medical device and complies with all the relevant Blood Regulations.



Simulated load and temperature probe



Hygienic, easy to read display and controller

Can be fitted with 21 stainless steel drawers designed to make storage, retrieval and inventory keeping easy without compromising air circulation.

Please order:

18 x BAPF2134A
3 x BAPF2134B



Part Number	RVPF2134MD
Capacity (Litres)	505
Style	Freestanding
Exterior Dimensions (HxWxD)	1880 x 770 x 722mm
Interior Dimensions (HxWxD)	1490 x 645 x 545mm
Weight	110kg
Energy Consumption	1.02 kWh/24hr
Operating Temperature	-34°C +/- 2°C
Interior Construction	Painted Steel
Exterior Construction	Powder Coated steel
Door Style	Solid
Door Hinging	Right
Open Door alarm	✓
Load Temperature Alarm	✓
Air Temperature Alarm	✓
Mains Failure Alarm	✓
Air Temperature Display	✓
Load Temperature Display	✓
7 Day Chart Recorder	✓
Lock	Manual Key
Battery Back up (Display)	✓
Shelves	6
Stainless steel drawers (optional extra)	21
Remote alarm contacts	✓
Access Port	✓
Auto Defrost	X
Internal Light	X
Eco Friendly Refrigerant	✓
Certified Medical Device in Accordance with Directive 93/42/EEC. Complies with BS4376:Part 2:1996	✓