

## **PRODUCT INFORMATION**

**Model: RLPR02042**

**Type: Labcold Sparkfree Laboratory Refrigerator**

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



**Capacity: 60 litres (gross)**

**Construction: White enamel exterior, fully moulded plastic interior**

**Shelves Supplied: 2, wire mesh**

**Exterior Dimensions (mm): 572 high x 485 wide x 520 deep**

**Features:**

**Sparkfree Interior**

**Door Lock**

**Easy-grip door handle**

**CFC and HFC free**

**Temperature set-point range: +2 to +10°C**

**Automatic Defrost**

**Order Reference: RLPR 02042**