

PRODUCT INFORMATION

Model: RLPR02042

Type: Labcold Sparkfree Laboratory Refrigerator

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



Capacity: 60 litres (gross)

Construction: White enamel exterior, fully moulded plastic interior

Shelves Supplied: 2, wire mesh

Exterior Dimensions (mm): 572 high x 485 wide x 520 deep

Features:

- Sparkfree Interior**
- Door Lock**
- Easy-grip door handle**
- CFC and HFC free**
- Temperature set-point range: +2 to +8°C**
- Automatic Defrost**

Order Reference: RLPR 02042