



Insulated Transport Box RTBB0006



A professional, low cost, reusable and wipe clean insulated box for transporting samples, specimens and biological materials.

Made from wipe clean plastics and coloured red, these boxes are supplied flat packed in quantities of 10 and come complete with carrying handles, adhesive labels and full assembly instructions.


Designed to keep contents cool while in transit between hospital departments until the contents can be placed in a suitable medical refrigeration unit. Each box has a thick white insulation layer to prevent the product from heating up.

Please note this product is not suitable for keeping blood products cool for extended periods of time.

- ❄ Low cost
- ❄ Fully insulated interior which can be removed
- ❄ Reusable
- ❄ Wipe clean exterior
- ❄ Light weight
- ❄ No cold packs needed
- ❄ Supplied in pack of 10
- ❄ Robust

Part Number	RTBB0006
Configuration	Portable
Capacity (bags)	6 (approx)
Height x Width x Depth (mm)	180 x 336 x 165
Weight (kg)	0.309
Pack Size	10

Easy self assembly



Flat Pack



Fold into box



Add insulated layer



Close the insulated layer



Add handle



Add label