



# LABCOLD™

medical & scientific refrigeration



## RLSD0600 Refrigerated Cooled Incubator

Part Number	RLSD0600
Capacity (Litres)	570
Style	Freestanding
Exterior Dimensions (HxWxD)	1980 x 720 x 840mm
Interior Dimensions (HxWxD)	1405 x 600 x 570mm
Weight	140kg
Internal Light	LED
Factory Set Point	+20°C
Operating Temperature	+5°C to +40°C
Interior Construction	Stainless Steel
Exterior Construction	Enamelled Steel
Door Style	Solid
Door Hinging	Right
Inner Doors	2
Shelves	4
Shelf dimension (WxDxH)	580 x 530 x 15mm
Spare shelves available	✓
Temperature display	✓
Microprocessor controlled temperature regulation	✓
Fan assisted temperature control	✓
Access Port	✓
Auto Defrost	X
Alarm	Optional Extra
Eco Friendly Refrigerant	✓

This 570 litre cooled incubator is ideal for microbiology or any application where precise and consistent temperatures at or below room temperature are required.

Digital Microprocessor Controller



With a robust stainless steel interior, enamelled steel exterior and fully adjustable shelves, this incubator is versatile and built to last.

Labcold Cooled Incubators are also available in 88 litre, 150 litre and 600 litre capacity

Most models have the choice of glass or solid door

01256 705 570

[www.labcold.com](http://www.labcold.com)

