



LABCOLD™

medical & scientific refrigeration



RLFR1804 Basic Laboratory Refrigerator

Part Number	RLFR1804
Capacity (Litres)	505
Style	Freestanding
Exterior Dimensions (HxWxD)	1710 x 765 x 725mm
Interior Dimensions (HxWxD)	1430 x 640 x 500mm
Temperature Display	Analogue
Energy Consumption	2.4 kWh/24h
Operating Temperature	+4°C
Temperature Set Point Range	+2°C to +10°C
Sparkfree	✗
Interior Construction	Painted Aluminium
Exterior Construction	White Enamel
Shelves	4
Door Hinging	Right
Door Style	Solid
Feet	2 feet at front 2 castors at back
Door Lock	✓
Auto Defrost	✓
Fan Circulated	✓
Ammonia Free	✓

Constructed from enamelled steel, this fridge from the **Labcold** Basic Range is designed to provide basic storage at a lower price while still providing the confidence of the **Labcold** name and support.

Lower cost production methods have been used in order to offer the best price possible without significant sacrifices in the important areas of function and reliability.



Analogue
Temperature
Display

This 505 litre fridge has four adjustable shelves and castors at the back of the unit to allow for easy manoeuvrability and ease of cleaning. A temperature display shows you the temperature inside the cabinet and has autodefrost for your convenience.

01256 705 570

www.labcold.com

