



Powered by IntelliCold™ Technology

The Labcold Breast Milk refrigerator has been specifically designed to store human milk in optimum conditions. With three adjustable shelves and an interior fan for even temperature distribution, this fridge is fitted with a lock for extra security.



The simple to use four button controller displays the internal temperature of the fridge, but it also shows minimum and maximum temperatures which can be displayed at the touch of a button. It also alarms if the internal temperature drops below 0°C or rises above 5°C. The Breast Milk fridge is also fitted with a door open alarm and, if batteries are fitted, a mains failure alarm.



The RLBM0510 comes with a self affix easy open handle designed to stand proud of surrounding units for easy opening and closing.

We are confident that this breast milk refrigerator will provide you with many years of fault free service. However, for extra peace of mind your RLBM0510 comes complete with a FREE 1 year parts and labour warranty!*

The controller automatically stores temperature data and fridge performance. A specifically encoded SD card is available as an optional extra for data transference so these records can be stored to your PC. Why not ask for details when you place your order?



RLBM0510 Breast Milk Refrigerator

Part Number	RLBM0510
Capacity (Litres)	150
Style	Under Bench
Exterior Dimensions (HxWxD)	820 x 600 x 610mm
Interior Dimensions (HxWxD)	690 x 510 x 410mm
Energy Consumption	0.66 kWh/24h
Operating Temperature	+3°C
Interior Construction	Painted Corrosion proof Aluminium
Exterior Construction	Powder Coated Steel
Shelves	3
Internal Drawers Available	✓
Door Style	Solid
Door Hinging	Right (can be reversed at time of manufacture)
Door Lock	✓
Digital Temperature Control	✓
Digital Temperature Display	✓
Internal Light	✓
High/Low Temperature alarm	✓
Minimum/Maximum Temperature Recording	✓
Door Open alarm	✓
Mains Failure alarm	✓
Alarm Battery Back-Up	✓
Off cycle auto defrost	✓
Fan Assisted Cooling	✓
Remote Alarm Contacts	✓
Eco Friendly Refrigerant	✓
Ammonia Free	✓

*Terms and conditions apply, please visit the website

01256 705 570

www.labcold.com

