



LABCOLD™

medical & scientific refrigeration



RLVF1514 Sparkfree Freezer

The Labcold Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this freezer.



For extra security all sparkfree models are fitted with a built in lock and digital temperature display. The display features high and low audio and visual temperature alarms as standard in addition to a door open alarm. The controller also records the minimum and maximum temperature since the last reset.

Part Number	RLVF1514
Capacity (Litres)	400
Style	Free standing
Exterior Dimensions (HxWxD)	1855 x 600 x 585mm
Interior Dimensions (HxWxD)	1620 x 470 x 480mm
Weight	79kg
Energy Consumption	1.9kWh/24hr
Operating Temperature	-12°C to -24°C
Factory Set Point	-20°C
Interior Construction	Moulded ABS
Exterior Construction	Painted Steel
Door Style	Solid
Door Hinging	Right (Can be reversed at time of manufacture)
Door Lock	✓
Door Open Alarm	✓
Fixed Shelves	6
Shelf Dimensions (W x D)	480 x 410mm
Access Port	Optional Extra
Feet	2 (adjustable)
Rollers	2 at rear
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Min/Max Recording	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓
2 Year Warranty (UK only)	✓

01256 705 570

www.labcold.com

