



RLCG0600 Refrigerated Cooled Incubator



This 570 litre cooled incubator is ideal for microbiology or any application where precise and consistent temperatures at or below room temperature are required.

Digital
Microprocessor
Controller



With a robust stainless steel interior, enamelled steel exterior and fully adjustable shelves, this incubator is versatile and built to last.

Part Number	RLCG0600
Capacity (Litres)	570
Style	Freestanding
Exterior Dimensions (HxWxD)	1990 x 710 x 800mm
Weight	170kg
Internal Light	LED
Temperature Uniformity	±0.5°C at 4°C ±1.0°C at 20°C
Operating Temperature	+5°C to +40°C
Interior Construction	Stainless Steel
Exterior Construction	Enamelled Steel
Door Style	Glass
Door Hinging	Right
Door Lock	✓
Inner Doors	2
Shelves	4
Shelf dimension (WxD)	500 x 620mm
Spare shelves available	✓
Temperature display	✓
Microprocessor controlled temperature regulation	✓
Fan assisted temperature control	✓
Access Port	✓
Auto Defrost	X
Alarm	Optional extra
Eco Friendly Refrigerant	✓

Labcold Cooled Incubators are also available in 150 litre and 340 litre capacity

Most models have the choice of glass or solid door