



RLVL03203 Basic Laboratory Freezer



Constructed from enamelled steel, this freezer from the **Labcold** Basic Range is designed to provide basic storage at a lower price while still providing the confidence of the **Labcold** brand and support.

Lower cost production methods have been used in order to offer the best price possible without significant sacrifices in the important areas of function and reliability.



One clear fronted flap.



Two clear plastic drawers.

Part Number	RLVL03203
Capacity (Litres)	70
Style	Underbench
Exterior Dimensions (HxWxD)	835 x 500 x 490mm
Interior Dimensions (HxWxD)	670 x 305 x 365mm
Energy Consumption	0.54 kWh/24h
Operating Temperature	-20°C
Temperature Set Point Range	-18°C to -23°C
Sparkfree	✗
Interior Construction	Moulded Plastic
Exterior Construction	Enamelled Steel
Drawers	2 + 1 flap
Door Hinging	Right
Door Style	Solid
Feet	2 adjustable feet at front
Manual Defrost	✓
Eco Friendly Refrigerant	✓
Ammonia Free	✓

This 70 litre freezer has two clear plastic drawers and one flap located at the top. Two adjustable feet ensure it remains stable and manual defrost provides you with greater operator control over your temperature sensitive contents.