



LABCOLD™

medical & scientific refrigeration



RPF05043 Solid Door Pharmacy & Vaccine Refrigerator*

Part Number	RPF05043
Capacity (Litres)	150
Style	Under Counter
Exterior Dimensions (HxWxD)	820 x 600 x 640 mm
Interior Dimensions (HxWxD)	660 x 545 x 355 mm
Energy Consumption	1.7 kWh/24h
Operating Temperature	+4°C
Temperature Set Point Range	+2°C to +8°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
Shelves	2
Shelf dimensions (W x D)	530 x 325 mm
Door Hinging	Right hand only
Door Style	Solid (Glass Door Model available RPF05043)
Internal Light	Glass Door Model only
Rear Castors	✓
Door Lock	✓
Access Port	Optional Extra
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature alarm	✓
Low Temperature alarm	✓
Door Open alarm	✓
Auto defrost	✓
Fan Assisted Cooling	✓
Environmentally friendly insulation	✓
Eco Friendly Refrigerant	✓
Ammonia Free	✓
DH Green Book, RPS/ GPhC, NaTHNaC and RCVS Guidelines compliant	✓

Robustly constructed from Stainless Steel this refrigerator from the **Labcold** Advanced range is designed to provide close temperature control for temperature sensitive products.



The fan circulated interior has no internal 'step', making it easy to clean and stock and to allow for maximum air flow for even temperature distribution. The two shelves are fully adjustable, so it is possible to maximise the space available.

Fitted with a sturdy door lock, this Advanced refrigerator is ideal for busy environments such as busy laboratories or pharmacies.

* Glass door model available. Please order RPF05043.

01256 705 570

www.labcold.com

