



Repeat Alarm RLMR0016

The **Labcold** Repeat Alarm can be remotely located up to 100 metres (wiring required) from the equipment it is monitoring. A self-recharging battery keeps you safe in the knowledge that the alarm will alert even in the event of a mains power failure providing you with crucial time to save your stock.



- ❄ Mains powered with self-recharging battery
- ❄ Mute button
- ❄ Red LED warning light
- ❄ Installed by **Labcold** Service Technicians (or suitably qualified person)

Part Number	RLMR0016
Configuration	Wall Mounted
Dimensions (H x W x D)	123 x 123 x 60mm
Connectivity	Twin Wire
Power Connections	240 volt AC, 1 amp (Fused spur recommended)