



**NEW Digital Lock
only for IntelliCold® Pharmacy Refrigerators
uniquely from**



Ergonomic sure touch buttons



*Sturdy metal lock engages when the door closes
(lock down for image purposes only)*

Fitting snugly to the front of the pharmacy fridge this NEW digital lock, unique to Labcold offers an even higher degree of safety and security for cold chain pharmaceutical storage.

The 10 digit combination lock can be set up by the user and the code only given to those who need it, meaning access is easily restricted. It can be changed when users change, allowing you to keep track of exactly who has access to a particular refrigerator.

Automatically locking when the door is shut, this lock also saves you time and prevents the risk of a pharmacy refrigerator being left unlocked.

Designed by Labcold to only fit on all our IntelliCold® Pharmacy Range, this sturdy digital lock, like our fridges, has been purpose designed to assist busy professionals store temperature sensitive medicines conveniently, securely and safely. Why risk your patients and pharmaceuticals when a new level of safety is at your fingertips?

**No more lost or misplaced keys
with the Labcold Digital Lock**

**You are in control over who has
access to the fridge**

**Eliminates the problem of
fridges being left unlocked**



**Safety at your
fingertips**

www.labcold.com

**Order DigilockDF10
with your IntelliCold® Pharmacy Fridge**