





Using tried and tested cascade refrigeration this new under bench ULT freezer from Labcold also features VIP insulation to keep the footprint even smaller.



Intuitive Control Panel

For sample security the LULT80100 has a door lock plus features all the relevant alarms. It also comes complete with factory fitted remote alarm contacts and an access port.

For extra peace of mind, this freezer has a FREE 5 year parts warranty and a 3 year labour warranty.



Labcold -80°C Ultra Low Temperature Freezer LULT80100

Part Number	LULT80100
Capacity (Litres)	100
Style	Under bench
Exterior Dimensions (HxWxD)	800 x 950 x 735 mm
Interior Dimensions (HxWxD)	500 x 450 x 450 mm
Weight	120 kg
Operating Temperature	-80°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
VIP insulation	\checkmark
Door Style	Solid
Door Hinging	Right
Open Door alarm	\checkmark
High Temperature Alarm	\checkmark
Mains Failure Alarm	\checkmark
Digital Temperature Display	\checkmark
Access Port	\checkmark
Remote Alarm Contacts	\checkmark
Key Operated Door Lock	\checkmark
Battery Back up (Display & Alarms only)	✓
Shelves	1
Inner Door	×
Rollers	\checkmark
Cascade Refrigeration	\checkmark
FREE 5 year parts and 3 year labour warranty (UK only)	\checkmark
Energy Consumption (kWh/24h)	8.7
Eco Friendly Refrigerant	✓

01256 705 570

www.labcold.com



