



LABCOLD™

medical & scientific refrigeration



Using tried and tested cascade refrigeration this new under bench ULT freezer from Labcold also features VIP insulation to keep the footprint even smaller.



Intuitive Control Panel

For sample security the LULT80100 has a door lock plus features all the relevant alarms. It also comes complete with factory fitted remote alarm contacts and an access port.

For extra peace of mind, this freezer has a FREE 5 year parts warranty and a 3 year labour warranty.



**Need racking?
Why not ask for a
quote?**

Labcold -80°C Ultra Low Temperature Freezer LULT80100

| Part Number | LULT80100 |
|--|--------------------|
| Capacity (Litres) | 100 |
| Style | Under bench |
| Exterior Dimensions (HxWxD) | 800 x 950 x 735 mm |
| Interior Dimensions (HxWxD) | 500 x 450 x 450 mm |
| Weight | 120 kg |
| Operating Temperature | -80°C |
| Interior Construction | Stainless Steel |
| Exterior Construction | Stainless Steel |
| VIP insulation | ✓ |
| Door Style | Solid |
| Door Hinging | Right |
| Open Door alarm | ✓ |
| High Temperature Alarm | ✓ |
| Mains Failure Alarm | ✓ |
| Digital Temperature Display | ✓ |
| Access Port | ✓ |
| Remote Alarm Contacts | ✓ |
| Key Operated Door Lock | ✓ |
| Battery Back up (Display & Alarms only) | ✓ |
| Shelves | 1 |
| Inner Door | ✗ |
| Rollers | ✓ |
| Cascade Refrigeration | ✓ |
| FREE 5 year parts and 3 year labour warranty (UK only) | ✓ |
| Energy Consumption (kWh/24h) | 8.7 |
| Eco Friendly Refrigerant | ✓ |

01256 705 570

www.labcold.com

