




Platelet Agitator and Racks

Labcold Platelet Agitators are manufactured in the UK and are certified as Medical Devices in accordance with UK MDR 2002. They have been engineered to agitate at the correct speed as defined by the World Health Organisation.


Labcold Platelet Agitators can accept either a six tray or an eight-tray rack. When used in conjunction with a Labcold Incubator the agitator stops when the door is opened and resumes movement when the door is closed.

To ensure the platelet bags do not come off the agitator racks, anti slip mats can be purchased.

All Labcold Agitators and Racks come complete with a full **FREE** 1 year parts and labour warranty.



Platelet Agitator - AGIT1006MD



Platelet Agitator + RACK1006



Platelet Agitator + RACK1008

Model	AGIT1006MD	RACK1006	RACK1008
Capacity (Depending on local pack size)	N/A	6	8
Height (mm)	145	197	253
Width (mm)	385	341	386
Depth (mm)	290	249	337
Weight (kg)	4.8	2.9	5.0
Exterior Construction	Painted Mild Steel	Stainless Steel/ Aluminium	Stainless Steel/ Aluminium
Colour	White	N/A	N/A
Style	Upright	Upright	Upright
Shelves	Rigid steel stand to support rack	6 slide out shelves	8 slide out shelves
Stroke Per Minute (Range)	40 - 60 rpm	N/A	N/A
Certified Medical Device in accordance with UK MDR 2002	✓	✓	✓