

RLPI0090MD Platelet Incubator (Medical Device)



Developed specifically for the safe storage of platelets at the optimum temperature. This Incubator is designed to ensure MHRA compliance.

Featuring a keypad lock for security and a built in chart recorder, this Incubator is packed with features that meet the BSQR standard for storing platelets.

Like the rest of the Labcold Blood Product Storage range this incubator is a certified medical device in accordance with EU Directive 93/42/EEC.

Part Number	RLPI0090MD
Capacity (Litres)	88
Style	Bench Top
Exterior Dimensions (HxWxD)	810 x 480 x 530mm
Interior Dimensions (HxWxD)	565 x 420 x 335mm
Energy Consumption	2.7 kWh/24hr
Operating Temperature	+22°C +/- 2°C
Interior Construction	Stainless Steel
Exterior Construction	Powder Coated steel
Door Style	Glass
Door Hinging	Right
Open Door alarm	✓
Load Temperature Alarm	✓
Air Temperature Alarm	✓
Mains Failure Alarm	✓
Motion Failure Alarm	✗
Air Temperature Display	✓
Load Temperature Display	✓
7 Day Chart Recorder	✓
Lock	Digital Keypad
Battery Back up (Display)	✓
Microprocessor Control	✓
Remote alarm contacts	✓
Access Port	✓
Auto Defrost	✗
Internal Light	✓
Eco Friendly Refrigerant	✓
Certified Medical Device in Accordance with Directive 93/42/EEC	✓



This deceptively spacious incubator has sufficient capacity to hold a Labcold flatbed Platelet Agitator (model number AGIT006MD) with a 6 tray rack which can be ordered separately.



Hygienic, easy to read display and controller



Simulated load alarm