

# Repeat Refrigerator Alarm **RLDR0001**

## *Intelligently Designed for Safety and Security*

Designed to work with the remote connections fitted as standard on the back of the Labcold IntelliCold® range of Pharmacy and Vaccine refrigerators, this battery operated repeat alarm allows you to hear and respond to a refrigerator alarm even if you are not in the same room. It is also ideal for noisy areas to boost the audible alarm fitted to the fridge, ensuring a rapid response.



Remote alarm contacts are located at the back of the fridge. The Repeat Refrigerator alarm kit is connected to the 2 pin connector on the left.

Supplied with a 5 meter long cable, which simply connects with the remote connections situated on the back of any refrigerator from the Labcold IntelliCold® range, the alarm is powered by a 9v battery, (not included). This means it works even if there is a power failure and no mains wiring is required, making it quick and easy to install.



Order **RLDR0001**

[www.labcold.com](http://www.labcold.com)

