



LABCOLD™

medical & scientific refrigeration



Labcold -50 to -86°C Ultra Low Temperature 390 litre Freezer

| Part Number | LULTO400U |
|---|---------------------|
| Capacity (Litres) | 390 |
| Style | Freestanding |
| Exterior Dimensions (HxWxD) | 1990 x 840 x 750mm |
| Interior Dimensions (HxWxD) | 1310 x 610 x 500mm |
| Weight | 280kg |
| Operating Temperature | -50°C to -86°C |
| Interior Construction | Stainless Steel |
| Exterior Construction | Epoxy Painted Steel |
| Door Style | Solid |
| Door Hinging | Right |
| Remote Alarm contacts | ✓ |
| High Temperature Alarm | ✓ |
| Low Temperature Alarm | ✓ |
| Mains Failure Alarm | ✓ |
| Door open alarm | ✓ |
| Low Battery Alarm | ✓ |
| Condenser High Temperature alarm | ✓ |
| Digital Temperature Display | ✓ |
| Key Operated Door Lock | ✗ |
| Password Protected Temperature Settings | ✓ |
| Battery Back up (Display and Alarms) | ✓ |
| Shelves | 3 |
| Inner Doors | 4 |
| Castors | 4 |
| Access Port | ✓ |
| Eco Friendly Refrigerant | ✓ |
| 5 year warranty 3 x years parts and labour 2 x years parts only | ✓ |

Labcold Ultra Low Temperature freezers use the tried and trusted cascade refrigeration and have been designed to occupy less than 1m² of floor space.



Complete with an easy to use digital controller conveniently situated on the front of the freezer, the temperature settings can be password protected and it always displays the set point for extra peace of mind. It also features alarm battery

back -up so should there be a mains power failure the freezer will alarm immediately to give you as much time as possible to deal with the emergency.

Backed by a 5 year parts and 3 year labour warranty, this Ultra Low Temperature freezer from Labcold is designed to be a robust and reliable place to store your samples.

We can also supply racking, why not ask for details at time of order?



01256 705 570

www.labcold.com

