

RLPR0517

Sparkfree Refrigerator



The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this fridge.



For extra security all sparkfree models are fitted with a built in lock and digital temperature display. The display features high and low audio and visual temperature alarms as standard in addition to a door open alarm.

Part Number	RLPR0517
Capacity (Litres)	150
Style	Under Bench
Exterior Dimensions (HxWxD)	835 x 600 x 600mm
Interior Dimensions (HxWxD)	690 x 480 x 410mm
Weight	45kg
Energy Consumption	1.1kWh/24hr
Operating Temperature	+0°C to +10°C
Factory Set Point	+5°C
Interior Construction	Powder Coated Steel
Exterior Construction	Powder Coated Steel
Door Style	Self Closing
Door Hinging	Right (can be reversed please ask at time of order)
Door Lock	✓
Shelves	3
Shelf dimension (WxD)	2 x 475 x 385mm 1 x 475 x 245mm
Spare shelves available	✓
Feet	4
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Door Open Alarm	✓
Access Port	✓
Off Cycle Auto Defrost	✓
Eco Friendly Refrigerant	✓
2 Year Warranty (UK only)	✓