



RLVF07203 Sparkfree Laboratory Freezer

Part Number	RLVF07203
Capacity (Litres)	170
Style	Freestanding
Exterior Dimensions (HxWxD)	1440 x 558 x 610mm
Interior Dimensions (HxWxD)	1240 x 390 x 340mm
Weight	51kg
Operating Temperature	-20°C
Interior Construction	Moulded Plastic
Exterior Construction	Powder Coated Steel
Sparkfree	✓
Shelves	1
Bins	4
Door Type	Solid
Door Hinging	Right
Defrost	Manual
High Temperature Alarm	Optional Extra
Low Temperature Alarm	Optional Extra
Feet	2
Rollers	2
Energy Consumption	0.801 kWh/24h
Eco Friendly Refrigerant	✓

Designed for the storage of potentially flammable contents, this sparkfree freezer from Labcold has all sources of ignition removed from the inside of the cabinet making it an important part of laboratory safety.

Constructed from powder coated steel with a moulded plastic interior, this freezer has been robustly constructed for professional use and is manual defrost only, which puts you in control of the temperature at all times.

Did you know that 25% of chemical accidents are caused by improper storage?

Source: US Chemical Safety and Laboratory Investigation Board 2010.

This freezer has four fixed shelves, including a top fixed shelf fitted with a plastic flap to prevent a large temperature variation when the door is opened.

Backed by the comprehensive Labcold 1 year parts and labour warranty, this freezer is made to order at our UK Headquarters and can be fitted with a temperature alarm as an optional extra.