



# LABCOLD™

medical & scientific refrigeration



## RLPR1517 Sparkfree Refrigerator

The Labcold Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

**Warning: Do not store open containers of volatile substances in this fridge.**



CLAIM  
YOUR  
**FREE**  
YEAR  
WARRANTY

For extra security all sparkfree models are fitted with a built in lock and digital temperature display. The display features high and low audio and visual temperature alarms as standard in addition to a door open alarm.

Part Number	RLPR1517
Capacity (Litres)	439
Style	Free standing
Exterior Dimensions (HxWxD)	1875 x 600 x 700mm
Interior Dimensions (HxWxD)	1724 x 500 x 535mm
Net Weight	93kg
Energy Consumption	2.55kWh/24hr
Operating Temperature	0°C to +10°C
Factory Set Point	+5°C
Interior Construction	Painted Aluminium
Exterior Construction	Painted Steel
Door Style	Solid
Door Hinging	Right (Can be reversed at time of manufacture)
Door Lock	✓
Door Open Alarm	✓
Adjustable Shelves	6
Shelf dimensions (WxD)	475mm x 470mm
Fan Circulated	✓
Rollers	4
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Access Port	✓
Automatic Defrost	✓
2 Year Warranty (UK only)	✓

01256 705 570

[www.labcold.com](http://www.labcold.com)

