




# LABCOLD™

medical & scientific refrigeration

## PHARMDRAWER1901 Divided Drawers for RLDF/G1019 and 1519



Designed exclusively to fit Labcold IntelliCold® pharmacy refrigerators 340 litres and over, these divided pharmacy drawers help keep your refrigerator neat and tidy enabling easy stock rotation.

Constructed from wipe clean plastic, these drawers can be retro fitted and can even be used in conjunction with shelves, giving you the choice on how best to utilize the space in your refrigerator. We



do, however, only recommend the use of 5 drawers in 340 litre models and 6 in 440 litre fridges.

Featuring an easy glide mechanism which ensures the product stays in place when the drawers are opened and closed, these drawers are designed to allow maximum air flow around the product.



Easy to move dividers



Use with shelves for maximum flexibility



Designed allow maximum air flow

01256 705 570

[www.labcold.com](http://www.labcold.com)

