

Master Key for RLDF/GDIGLOCK IntelliCold[®] Pharmacy Refrigerators

LABLOCK3KEYS



Labcold IntelliCold[®] Digital Lock closed



Labcold IntelliCold[®] Digital Lock open

All Labcold IntelliCold[®] Pharmacy Refrigerators now come with the option of a digital lock. These are the models with DIGLOCK in the product code.

Digital locks solve the problem of lost keys and generic pharmacy refrigerator keys so it is easier to control who has access to the refrigerator. Plus, if a member of the team leaves, you can simply change the code to maintain security.

The digital lock is hardwired in to the mains electric of the refrigerator, so you never have to worry about changing a battery. Batteries in the Labcold IntelliCold[®] fridges enable not only the mains fail alarm and continue to record temperature data in the event of a power failure but they generate enough power to allow you to open the digital lock should you experience a power cut.

However, what happens if the fridge is 'retrieved' or repurposed? Sometimes the code isn't passed on to those who need it. That's why the new digital lock has a key lock and you can purchase a master key to ensure you are never locked out of our refrigerators. Especially useful for Estates Department, the keys are made for Labcold locks only for extra security.

**Ensures that the digital lock
can be opened, even if no one
remembers the code!**

Ideal for Estates Departments

