



# LABCOLD™

medical & scientific refrigeration



## RAFR05204 Advanced Under Counter Freezer

Part Number	RAFR05204
Capacity (Litres)	115
Style	Under Counter
Exterior Dimensions (HxWxD)	835 x 600 x 600mm
Energy Consumption	1.5kWh/24h
Operating Temperature	-20°C
Temperature Set Point Range	-18°C to -22°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
Shelves	2
Shelf dimensions (W x D)	460 x 4075mm 460 x 240mm
Door Hinging	Right (can be reversed at time of order)
Self Closing Door	✓
Door Lock	✓
Levelling Feet	4
Access Port	✓
Remote Alarm Contacts	✓
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature alarm	✓
Low Temperature alarm	✓
Door Open alarm	✓
Manual defrost	✓
GWP Friendly Refrigerant	✓
Ammonia Free	✓
2 Year Warranty*	✓

Robustly constructed with a stainless Steel exterior this freezer from the **Labcold** Advanced range is designed to provide close temperature control for temperature sensitive products .



Insulated with Polyurethane Water Base Insulation the refrigerant is also ammonia, HCFC and CFC free and uses a GWP friendly refrigerant.

Fitted with a sturdy door lock, this Advanced freezer is ideal for demanding environments such as busy laboratories and features a manual defrost putting you in control.

Remote alarm contacts and an access port are also fitted as standard.

\*Parts only warranty outside of the UK

01256 705 570

www.labcold.com

