



LABCOLD™

medical & scientific refrigeration



RLPL04044 Basic Laboratory Refrigerator

| Part Number | RLPL04044 |
|-------------------------------------|--|
| Capacity (Litres) | 92 |
| Style | Underbench |
| Exterior Dimensions (HxWxD) | 842 x 475 x 448mm |
| Weight | 22kg |
| Energy Consumption | 0.29 kWh/24h |
| Noise Level | 39 dB(A) |
| Operating Temperature | +4°C |
| Temperature Set Point Range | +2°C to +10°C |
| Interior Construction | Moulded Plastic |
| Exterior Construction | Enamelled Steel |
| Shelves | 2 + 1 bin |
| Sparkfree | ✗ |
| Door Hinging | Right (can be reversed at time of ordering) |
| Door Style | Solid |
| Feet | 2 adjustable feet at front |
| Internal Light | ✓ |
| Auto Defrost | ✓ |
| GWP Friendly Refrigerant | ✓ |
| Ammonia Free | ✓ |
| 1 year UK Parts and Labour Warranty | ✓ |

Constructed from enamelled steel, this refrigerator from the **Labcold** Basic Range is designed to provide basic storage at a lower price while still providing the confidence of the **Labcold** brand and support.

The fridge has two non adjustable shelves and a polycarbonate drawer at the bottom of the refrigerator. The adjustable feet ensures the fridge can be stabilised and with auto defrost and an internal light, the refrigerator has some useful convenience features.

To put you in control, the temperature can be easily adjusted with a dial located on the right hand side of the interior. Labcold will by default set the set point to +4°C.

Please note, this refrigerator is not sparkfree and can not have a lock or digital temperature display fitted to it.

This Basic Refrigerator is ideal for economic cold storage of items that are not volatile and where precise storage temperature is not product critical. As you would expect from Labcold, this Basic Refrigerator comes complete with a full one years parts and labour warranty.*



* Labour UK only