



# LABCOLD™

medical & scientific refrigeration

## RLVL03204 Basic Laboratory Freezer



Constructed from enamelled steel, this freezer from the **Labcold** Basic Range is designed to provide basic storage at a lower price while still providing the confidence of the **Labcold** brand and support.

The freezer has two clear plastic drawers, one which is the interior width of the freezer and a half drawer at the bottom. There is a storage compartment that runs the whole depth of the freezer with a spring loaded flap to keep contents safe at the top. The adjustable feet ensures the freezer can be stabilised and manual defrost provides you with greater operator control over your temperature sensitive contents.

To put you in control, the temperature can be adjusted with an easy dial located on the right hand side of the interior, however, Labcold will by default set the set point to -20°C.

Please note, this freezer is not sparkfree and can not have a lock or digital temperature display fitted to it.

This Basic Freezer is ideal for economic cold storage of items that are not volatile and where precise storage temperature is not product critical. As you would expect from Labcold, this freezer comes backed with a full one years parts and labour warranty.\*



Part Number	RLVL03204
Capacity (Litres)	61
Style	Underbench
Exterior Dimensions (HxWxD)	849 x 475 x 448mm
Energy Consumption	0.54 kWh/24h
Noise level	39dB(A)
Temperature Control	✓
Operating Temperature	-20°C
Temperature Set Point Range	-18°C to -23°C
Sparkfree	✗
Interior Construction	Moulded Plastic
Exterior Construction	Enamelled Steel
Drawers	2 + 1 flap
Door Hinging	Right (can be reversed at time of ordering)
Door Style	Solid
Feet	4 adjustable
Manual Defrost	✓
GWP Friendly Refrigerant	✓
Ammonia Free	✓
1 year UK Parts and Labour Warranty	✓

\* Labour UK only

01256 705 570

[www.labcold.com](http://www.labcold.com)

