



# LABCOLD™

medical & scientific refrigeration



## RPF05044 Advanced Pharmacy Under Counter Refrigerator

Part Number	RPF05044
Capacity (Litres)	150
Style	Under Counter
Exterior Dimensions (HxWxD)	835 x 600 x 600 mm
Energy Consumption	1.8 kWh/24h
Operating Temperature	+5°C
Temperature Set Point Range	+2°C to +8°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
Shelves	2
Shelf dimensions (W xD)	475 x 385mm 475 x 245mm
Door Hinging	Right hand (can be reversed at time of ordering)
Door style	Solid
Internal Light	✓
Door Lock	✓
Access Port	✓
Remote Alarm Contacts	✓
Digital Temperature Display	✓
High Temperature alarm	✓
Low Temperature alarm	✓
Door Open alarm	✓
Auto defrost	✓
Fan Assisted Cooling	✓
Eco Friendly Refrigerant	✓
FREE manufacturers' calibration certificate from our UKAS temperature calibration laboratory #8898	✓
2 Year UK Warranty	✓

Robustly constructed from stainless steel this refrigerator from the **Labcold** Advanced pharmacy range is designed to provide close temperature control for temperature sensitive vaccines, medicines and reagents.



The fan circulated interior ensures that there are no hot or cold spots and the robust construction means it's ideal for busy pharmacies where the door is frequently opened. Fitted with a sturdy door lock for extra security, remote alarm contacts and an access port are also fitted as standard.

For extra peace of mind, this Advanced Pharmacy refrigerator comes complete with a calibration certificate from our UKAS temperature calibration laboratory #8898 and is backed by a 2 year parts and labour warranty.\*

\*Parts only warranty outside of the UK

01256 705 570

[www.labcold.com](http://www.labcold.com)

