



Inspection and Calibration contracts

Labcold Inspection and Calibration contracts are designed to ensure single medical and laboratory refrigeration owners can get the same expert service as owners of our certified Medical Device products.


All pharmacy refrigerators that store vaccines must be calibrated annually and have a maintenance contract (DH Green Book, Chapter 3v2_1), so we designed Labcold's budget friendly inspection and calibration contract to do just this.

Our FGas trained engineer will carry out a full inspection of your fridge. If it's one of ours, because our engineers carry OEM spare parts routinely, any issues can be resolved quickly. They will also calibrate the temperature display, and because Labcold is an ISO 17025 accredited UKAS temperature calibration laboratory #8898 you can be sure that the procedure and the equipment used will provide a result that is repeatable and verifiable to international standards.



Typical Inspection and Calibration Checks

- ❄️ All settings checked and adjusted if incorrect
- ❄️ All control functions tested
- ❄️ Thorough inspection for any visible damage
- ❄️ Test to ensure all functions are working correctly
- ❄️ Calibration will be performed to prove that the temperature shown on the display is accurate
- ❄️ A calibration certificate will be issued by our UKAS accredited temperature calibration laboratory #8898



Our engineers covered the whole of the UK, including the Scottish and Channel Islands. Being a UKAS temperature calibration laboratory, you can be confident that they can calibrate data loggers and pharmacy refrigerators from any manufacturer to the same exacting standards. Because we are the UK's largest manufacturer of medical and scientific refrigeration, our friendly engineers will be

able to answer any queries you have about your fridge and be able to give you tips on how to run it better for longer. After all if you've seen them being made you know exactly how they work!

To find our more, contact our sales team at: contracts@labcold.com

