



## RLPR0217- ATEX Bench Top Sparkfree ATEX Refrigerator

The Labcold ATEX refrigerator range has been designed with an interior which is licensed as a Zone 2 hazardous area under 2014/34/EU (ATEX Directive) in order to store products which are classed in explosion group IIB+H2 and temperature class T6 under 2014/34/EU (ATEX Directive).

**Warning: Do not store open containers of volatile substances in this fridge.**



CLAIM  
YOUR  
**FREE**  
YEAR  
WARRANTY

All Labcold sparkfree ATEX refrigerators are fitted with a cylinder lock with two keys and a controller with a digital temperature display. The controller also features a high/low audio and visual temperature alarm in addition to a door open alarm.

Part Number	RLPR0217-ATEX
Capacity (Litres)	66
Style	Bench Top
Exterior Dimensions (HxWxD)	740 x 450 x 510mm
Interior Dimensions (HxWxD)	590 x 333 x 300mm
Weight	34kg
Energy Consumption	1.4kWh/24hr
Operating Temperature	+0°C to +10°C
Factory Set Point	+5°C
Interior Construction	Painted Aluminium
Exterior Construction	Powder Coated Steel
Door Style	Self Closing
Door Hinging	Right (can be reversed please ask at time of order)
Door Lock	✓
Door Earth Connection	✓
Condensate Collection	✓
Shelves	3
Shelf dimension (WxD)	2 x 346 x 203mm 1 x 346 x 150mm
Spare shelves available	✓
Feet	4
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Door Open Alarm	✓
Off cycle auto Defrost	✓
Low GWP Refrigerant	✓
2 Year Warranty (UK only)	✓