



## RLPR1517-ATEX Free Standing Sparkfree ATEX Refrigerator

The **Labcold** ATEX refrigerator range has been designed with an interior which is licensed as a Zone 2 hazardous area under 2014/34/EU (ATEX Directive) in order to store products which are classed in explosion group IIB+H2 and temperature class T6 under 2014/34/EU (ATEX Directive).

**Warning: Do not store open containers of volatile substances in this fridge.**



CLAIM  
YOUR  
**FREE**  
YEAR  
WARRANTY

All Labcold sparkfree ATEX refrigerators are fitted with a cylinder lock with two keys and a controller with a digital temperature display. The controller also features a high/low audio and visual temperature alarm in addition to a door open alarm.

Part Number	RLPR1517-ATEX
Capacity (Litres)	439
Style	Free standing
Exterior Dimensions (HxWxD)	1875 x 600 x 700mm
Interior Dimensions (HxWxD)	1724 x 500 x 535mm
Net Weight	93kg
Energy Consumption	2.55kWh/24hr
Operating Temperature	0°C to +10°C
Factory Set Point	+5°C
Interior Construction	Painted Aluminium
Exterior Construction	Painted Steel
Door Style	Solid
Door Hinging	Right (Can be reversed at time of manufacture)
Door Lock	✓
Door Earth Connection	✓
Condensate Collection	✓
Door Open Alarm	✓
Adjustable Shelves	6
Shelf dimensions (WxD)	475mm x 470mm
Fan Circulated (ATEX fan)	✓
Rollers	4
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Automatic Defrost	✓
Low GWP Refrigerant	✓
2 Year Warranty (UK only)	✓