



LABCOLD™

medical & scientific refrigeration



RLVF0417-ATEX Sparkfree ATEX Under Counter Freezer

Part Number	RLVF0417 - ATEX
Capacity (Litres)	124
Style	Under Bench
Exterior Dimensions (HxWxD)	835 x 600 x 600mm
Interior Dimensions (HxWxD)	680 x 460 x 410mm
Weight	45kg
Energy Consumption	0.8kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Powder Coated Steel
Exterior Construction	Powder Coated Steel
Door Style	Self Closing
Door Hinging	Right (can be reversed please ask at time of order)
Shelves	3
Shelf Dimensions (w x d)	2 of 460 x 407mm 1 of 460 x 240mm
Door Lock	✓
Door Open Alarm	✓
Door Earth Connection	✓
Feet	4
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓
2 Year Warranty (UK only)	✓

The Labcold ATEX freezer range has been designed with an interior which is licensed as a Zone 2 hazardous area under 2014/34/EU (ATEX Directive) in order to store products which are classed in explosion group IIB+H2 and temperature class T6 under 2014/34/EU (ATEX Directive).

Warning: Do not store open containers of volatile substances in this freezer.



CLAIM
YOUR
FREE
YEAR
WARRANTY

All Labcold sparkfree ATEX freezers are fitted with a cylinder lock with two keys and a controller with a digital temperature display. The controller also features a high/low audio and visual temperature alarm in addition to a door open alarm.

01256 705 570

www.labcold.com

