



# LABCOLD™

medical & scientific refrigeration



## RLVF1517-ATEX Free Standing Sparkfree ATEX Freezer

Part Number	RLVF1517- ATEX
Capacity (Litres)	406
Style	Free standing
Exterior Dimensions (HxWxD)	1875 x 600 x 700mm
Interior Dimensions (HxWxD)	1690 x 440 x 500mm
Net Weight	97kg
Energy Consumption	3.16kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Painted Aluminium
Exterior Construction	Painted Steel
Door Style	Solid
Door Hinging	Right
Door Earth Connection	✓
Door Lock	✓
Door Open Alarm	✓
Adjustable Shelves	6
Shelf dimensions (WxD)	413mm x 438mm
Fan Circulated (ATEX fan)	✓
Rollers	4
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Automatic Defrost	✗
2 Year Warranty (UK only)	✓

The Labcold ATEX freezer range has been designed with an interior which is licensed as a Zone hazardous area under 2014/34/EU (ATEX Directive) in order to store products which are classed in explosion group IIB+H2 and temperature class T6 under 2014/34/EU (ATEX Directive).

**Warning: Do not store open containers of volatile substances in this freezer.**

CLAIM  
YOUR  
**FREE**  
YEAR  
WARRANTY



All Labcold sparkfree ATEX freezers are fitted with a cylinder lock with two keys and a controller with a digital temperature display. The controller also features a high/low audio and visual temperature alarm in addition to a door open alarm.

01256 705 570

[www.labcold.com](http://www.labcold.com)

